

Gas Required for Spectrometers



Overview

Gas Supply Requirements Ensure that the laboratory site meets the gas requirements as listed in the following table. NOTE If the location of your gas supply.

Key Takeaways

Mass spectrometers require high-purity gases, typically nitrogen or helium, delivered at specific pressures and free from contaminants to ensure accurate and reliable operation.

- **Types of Gases**
 - Nitrogen (N₂): Commonly used as a source gas, collision gas, and C-Trap bath gas in instruments like the Orbitrap Astral. Ultra-high purity (UHP, 99.999%) nitrogen is mandatory for optimal performance .
 - Helium (He): Often used as a carrier gas in gas chromatography-mass spectrometry (GC-MS) systems for sample transport and separation .
 - Other gases: Depending on the instrument, high-purity argon or hydrogen may be used for specific ionization or collision processes .

Purity Requirements

- Ultra-high purity (UHP, 99.999%) is required for critical functions such as reagent carrier gas and collision gas .
 - High purity (HP) gases may be acceptable for less critical functions, but contaminants like particulates, oils, or moisture can damage the instrument .
 - Avoid using particulate filters in gas lines, as they can introduce contamination .
- ### Pressure and Delivery
- Typical nitrogen pressure for Orbitrap instruments: 600 ± 50 kPa (87 ± 7 psi), with a maximum line pressure of 650 kPa (94 psi), .

- Gas lines must be clean, free of particulates and oils, and made of suitable materials such as copper or stainless steel for UHP gases .
 - If the laboratory gas supply serves multiple instruments, a pre-regulator should be installed to maintain stable pressure for the spectrometer .
 - Tubing length and connection fittings (e.g., 1/8" Swagelok) should be compatible with the instrument and allow easy access for installation and maintenance .
- Consumption Considerations
- Gas consumption varies depending on the type of analysis and instrument usage. Continuous operation can require significant volumes of gas .

- Ensure the supply system can maintain the required pressure and purity even when other laboratory equipment is using the same gas source .

Summary

For reliable spectrometer operation, laboratories must provide high-purity gases at controlled pressures, with clean, compatible delivery lines. Nitrogen and helium are the most common gases, with UHP nitrogen being essential for collision and carrier functions. Proper installation, monitoring, and maintenance of gas lines are critical to prevent contamination and ensure accurate mass spectrometry results .

Article Content

Gas Supply Requirements

Ensure that the laboratory site meets the gas requirements as listed in the following table. NOTE If the location of your gas supply

Gas Requirements

What should you consider when selecting gases for your new GC or GC/MS system? Which gases will you be using? Did you c

Nitrogen gas

The mass spectrometer requires nitrogen gas for its different components. The required gas pressure is 600 ± 50 kPa (87 ± 7 psi),

Agilent AA Spectrometer Systems Site Preparation Checklist

The VGA gas outlet should be fitted with a shutoff valve to prevent the loss of approximately 50 mm per minute of inert gas, which

Gas Analysis by Mass Spectrometry

Small gas volume and high precision requirements (quantitation and m/z resolution) make dual-focusing sector systems the best

20.3: Mass Spectrometers

As we noted in Chapter 11, a mass spectrometer must operate under a vacuum to ensure that ions can travel long distances without

4.3: Mass Spectrometers

Figure 4.3.1 provides a general view of a mass spectrometer, and the requirements are described briefly below. Note that the

Gas Supply to Your Mass Spectrometer: Focus on Tubing

Choosing the correct tubing to supply gas to your mass spectrometer is an important decision. Discover what you should consider

A short guide to gas monitoring in LABORATORIES

Bear in mind that it is not a legal requirement to follow standards, but even if it is not a legal requirement to have a gas monitor; it is

Providing Gases for LC-MS Systems in a Safe, Convenient, and Cost ...

A liquid chromatographic system with a mass spectrometric detector (LC-MS) requires a supply of nitrogen as the

Spectrometer

A spectrometer is any instrument used to view and analyze a range (or a spectrum) of a given characteristic for a substance (e.g., a

Standard Operating Procedure (SOP)

Standard Operating Procedure (SOP) This Standard Operating Procedure (SOP) describes basic chemical safety information for

Introduction to mass analyzers : Shimadzu Scientific Instruments

Mass spectrometers are now used for an extremely diverse range of applications, each with its own characteristics. Therefore, it is

Gas Source Mass Spectrometry: Stable Isotope Geochemistry

Gas source mass spectrometers are the principal instruments used for measuring isotopic ratios of light elements, including

What is Mass Spectrometry? | Broad Institute

Mass spectrometry is an analytical tool useful for measuring the mass-to-charge ratio (m/z) of one or more molecules present in a

Mass Spectrometry Tutorial (Dr. Kamel Harrata) | Chemical ...

This tutorial discusses basic aspects of mass spectrometry that will be helpful to you in deciding the proper techniques and

Mass spectrometry

Mass spectrometry (MS) is an analytical technique that is used to measure the mass-to-charge ratio of

Spectrometers for Elemental Spectrochemical Analysis, Part I: The

Spectrometers for Elemental Spectrochemical Analysis, Part I: The Basic Spectrometer An overview of the

QTRAP Compendium

QTRAP® Technology Explained What makes a QTRAP system so unique? While having the capability to function like a standard

Spectrometers – Visual Encyclopedia of Chemical Engineering

Mass spectrometers are widely used in laboratories for educational purposes and in the field to study the emission of molecules.

4.11: Mass Spectrometry

Typically, liquid chromatography systems are paired with ESI-quadrupole mass spectrometers to take advantage of

Gas chromatography mass spectrometry basic principles | Agilent

Learn about the basic principle of gas chromatography mass spectrometry (GC/MS). This GCMS primer explains the fundamentals

Nitrogen gas and Liquid nitrogen for your NMR instrument

In order to work, NMR spectrometers require nitrogen or liquid nitrogen to maintain the temperature of the superconducting magnet.

How the Mass Spectrometer Works

How a mass spectrometer works An outline of what happens in a mass spectrometer Understanding what's going on The need for a

Impact of active gas purge on FTIR measurements

The results unequivocally demonstrate that gas purge is a simple yet highly effective and efficient method to remove water vapor and

General Chapters: MASS SPECTROMETRY

Gas and liquid chromatographs are widely used as sample inlet devices for mass spectrometers. These chromatographs provide for

Atomic Absorption Spectroscopy Overview

Atomic absorption spectroscopy is an analytical technique used for trace metals analysis. What is atomic absorption spectroscopy

A Comprehensive Guide to Mass Spectrometers: How They Work and ...

A mass spectrometer is a key tool for analyzing many substances. It has an ion source, a mass analyzer, and a detector. These

Mass spectrometry | Definition, Applications, Principle, & Facts ...

Mass spectrometry, analytic technique by which chemical substances are identified by the sorting of gaseous ions in

Gas Requirements

Gas requirements are specified in Site Preparation Checklists and Guides. This information is also included in the free Agilent

Spectrometer

A spectrometer is any instrument used to view and analyze a range (or a spectrum) of a given characteristic for a substance (e.g., a

Standard Operating Procedure (SOP)

Standard Operating Procedure (SOP) This Standard Operating Procedure (SOP) describes basic chemical safety information for

Introduction to mass analyzers : Shimadzu Scientific Instruments

Mass spectrometers are now used for an extremely diverse range of applications, each with its own characteristics. Therefore, it is

Gas Source Mass Spectrometry: Stable Isotope Geochemistry

Gas source mass spectrometers are the principal instruments used for measuring isotopic ratios of light elements, including

What is Mass Spectrometry? | Broad Institute

Mass spectrometry is an analytical tool useful for measuring the mass-to-charge ratio (m/z) of one or more molecules present in a

Mass Spectrometry Tutorial (Dr. Kamel Harrata) | Chemical ...

This tutorial discusses basic aspects of mass spectrometry that will be helpful to you in deciding the proper techniques and

The PinAAcle 900 Spectrometer Series : Quote, RFQ,

The PinAAcle™ 900 atomic absorption (AA) series, engineered with a range of innovative technological

Mass Spectrometry Basics

Mass spectrometers are often combined with equipment that separates compounds in mixtures before the compounds are

Gas Chromatography

Gas chromatography is a term used to describe the group of analytical separation techniques used to analyze volatile substances in

Providing Gases for LC-MS Systems in a Safe, Convenient, and Cost ...

In recent years, the various components of gas generation systems have been integrated into a single system to

Mass spectrometry

Mass spectrometry (MS) is an analytical technique that is used to measure the mass-to-charge ratio of ions. The results are

Spectrometers for Elemental Spectrochemical Analysis, Part I: The

Spectrometers for Elemental Spectrochemical Analysis, Part I: The Basic Spectrometer An overview of the

QTRAP Compendium

QTRAP® Technology Explained What makes a QTRAP system so unique? While having the capability to function like a standard

Spectrometers – Visual Encyclopedia of Chemical Engineering

Mass spectrometers are widely used in laboratories for educational purposes and in the field to study the emission of molecules.

Gas Supply Requirements

Ensure that the laboratory site meets the gas requirements as listed in the following table. NOTE If the location of your gas supply

Gas Requirements

What should you consider when selecting gases for your new GC or GC/MS system? Which gases will you be using? Did you c

Nitrogen gas

The mass spectrometer requires nitrogen gas for its different components. The required gas pressure is 600 ± 50 kPa (87 ± 7 psi),

Agilent AA Spectrometer Systems Site Preparation Checklist

The VGA gas outlet should be fitted with a shutoff valve to prevent the loss of approximately 50 mm per minute of inert gas, which

Gas Analysis by Mass Spectrometry

Small gas volume and high precision requirements (quantitation and m/z resolution) make dual-focusing sector systems the best

20.3: Mass Spectrometers

As we noted in Chapter 11, a mass spectrometer must operate under a vacuum to ensure that ions can travel long distances without

4.3: Mass Spectrometers

Figure 4.3.1 provides a general view of a mass spectrometer, and the requirements are described briefly below. Note that the

Gas Supply to Your Mass Spectrometer: Focus on Tubing

Choosing the correct tubing to supply gas to your mass spectrometer is an important decision. Discover what you should consider

A short guide to gas monitoring in LABORATORIES

Bear in mind that it is not a legal requirement to follow standards, but even if it is not a legal requirement to have a gas monitor; it is

Providing Gases for LC-MS Systems in a Safe, Convenient, and Cost ...

A liquid chromatographic system with a mass spectrometric detector (LC-MS) requires a supply of nitrogen as the

Gas Supply Requirements

Ensure that the laboratory site meets the gas requirements as listed in the following table. NOTE If the location of your gas supply

Gas Requirements

What should you consider when selecting gases for your new GC or GC/MS system? Which gases will you be using? Did you c

Nitrogen gas

The mass spectrometer requires nitrogen gas for its different components. The required gas pressure is 600 ± 50 kPa (87 ± 7 psi),

Agilent AA Spectrometer Systems Site Preparation Checklist

The VGA gas outlet should be fitted with a shutoff valve to prevent the loss of approximately 50 mm per minute of inert gas, which

Gas Analysis by Mass Spectrometry

Small gas volume and high precision requirements (quantitation and m/z resolution) make dual-focusing sector systems the best

20.3: Mass Spectrometers

As we noted in Chapter 11, a mass spectrometer must operate under a vacuum to ensure that ions can travel long distances without

4.3: Mass Spectrometers

Figure 4.3.1 provides a general view of a mass spectrometer, and the requirements are described briefly below. Note that the

Gas Supply to Your Mass Spectrometer: Focus on Tubing

Choosing the correct tubing to supply gas to your mass spectrometer is an important decision. Discover what you should consider

A short guide to gas monitoring in LABORATORIES

Bear in mind that it is not a legal requirement to follow standards, but even if it is not a legal requirement to have a gas monitor; it is

Providing Gases for LC-MS Systems in a Safe, Convenient, and Cost ...

A liquid chromatographic system with a mass spectrometric detector (LC-MS) requires a supply of nitrogen as the

Contact Us

Continue your research online:

Website: <https://alfa-pro.co.za>

Technical inquiry: <https://alfa-pro.co.za/contact.html>

This document is for preliminary research. Verify specifications against current standards, product documentation and project requirements.

